- (en) Charge completely the battery of the Meter before its first use.
- (es) Cargue la batería del Monitor por completo antes del primer uso.
- (de) Laden Sie die Batterie des Messgeräts vor dem ersten Gebrauch vollständig auf.
- (nl) De batterij van de meter moet vóór het eerste gebruik volledig worden opgeladen.
- (it) Caricare completamente la batteria del dispositivo prima del suo primo utilizzo.
- (fr) Charger complètement la batterie de l'appareil avant la première utilisation.
- (pt) Carregue totalmente a bateria do Monitor antes da sua primeira utilização.
- (cs) Před prvním použitím měřicího přístroje baterii zcela nabijte.
- (sl) Pred prvo uporabo baterijo Merilnika napolnite v celoti.

iLine Microsystems S.L.

Paseo Mikeletegi, 69
20009 Donostia - Gipuzkoa (SPAIN)

www.ilinemicsystems.com

© 2020 iLine Microsystems. microINR is trademark of iLine Microsystems.
The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by iLine Microsystems is under license.

0123

IND0001AA - Rev 2020-05

**friendly
technology**

microINR[®] link

- (en) Instructions for use
- (es) Instrucciones de uso
- (de) Bedienungsanleitung
- (nl) Gebruiksaanwijzing
- (it) Istruzioni per l'uso
- (fr) Instructions d'utilisation
- (pt) Instruções de utilização
- (cs) Návod k použití
- (sl) Navodila za uporabo

en

No part of these instructions for use or the products it describes may be reproduced by any means or in any form without prior consent in writing from iLine Microsystems.

The information contained in these instructions for use was correct at the time of printing. However, iLine Microsystems reserves the right to introduce changes to the specifications, equipment and maintenance procedures of the system at any time, without prior notification. Any substantial changes will be added to the next revision of these instructions for use. In case of conflict between the information of the instructions for use of the microINR Link Meter and the microINR Chip, the information contained in the instructions for use of the microINR Chip will have priority.

microINR is trademark of iLine Microsystems.

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by iLine Microsystems is under license.

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For monitoring of oral anticoagulation treatment
based on Vitamin K antagonist drugs.

CONTENTS	
1. INTRODUCTION	2
1.1 INTENDED USE	2
1.2 BEFORE STARTING TO USE THE microINR® SYSTEM	2
Safety information	2
1.3 ORAL ANTICOAGULANT THERAPY	4
Prothrombin Time and INR	4
1.4 MEASURING PRINCIPLE	5
2. microINR® SYSTEM	6
2.1 DESCRIPTION OF THE microINR® Link KIT	6
2.2 PARTS OF THE microINR® Link METER	6
2.3 CHARGING THE microINR® Link METER	6
2.4 SETTING THE TIME & DATE	7
2.5 PATIENT IDENTIFICATION (OPTIONAL)	7
2.6 QUALITY CONTROL	8
3. CONDUCTING THE TEST	8
3.1 PREPARING THE NECESSARY MATERIALS	8
3.2 MEASUREMENT PROCEDURE	9
Turning the Meter On	9
Inserting the Chip	9
Modes of Use	9
Conducting the Test	10

Test Result and Assay End	10
Turning the Meter Off	11
3.3 OBTAINING AND APPLYING THE CAPILLARY BLOOD SAMPLE	11
3.4 INTERPRETING THE RESULTS	12
3.5 LIMITATIONS OF USE	13
4. MEMORY	13
5. BLUETOOTH® FUNCTIONALITY	14
5.1 BLUETOOTH® CONFIGURABLE FEATURES	14
5.2 microINR® Link METER BLUETOOTH® MENU	14
5.3 DEVICE PAIRING	15
5.4 DEVICE UNPAIRING	16
5.5 TRANSFERRING DATA	17
5.6 SENDING ALL RESULTS	17
6. CLEANING AND DISINFECTING THE microINR® Link METER	18
7. ERROR GUIDE	19
8. ADDITIONAL INFORMATION	20
8.1 SPECIFICATIONS	20
8.2 WARRANTY	20
8.3 TECHNICAL SERVICE	20
8.4 SYMBOLS	21
8.5 GLOSSARY	23

1. INTRODUCTION

1.1 INTENDED USE

The microINR System is intended to monitor oral anticoagulation therapy (OAT) with vitamin K antagonist drugs. The microINR System determines quantitative prothrombin time (PT) in INR (International Normalized Ratio) units with fresh capillary blood performed by fingersticking.

The microINR System is a medical device for *in-vitro* diagnostics intended for professional and self-testing use.

1.2 BEFORE STARTING TO USE THE microINR® SYSTEM

These instructions for use will guide you, with comprehensive information, on the handling and use of the microINR Link Meter. Please, read them as well as the instructions for use of the microINR Chip carefully. Furthermore, do not forget to read the instructions for use of the disposable lancets and/or lancing device used to obtain the capillary blood sample.

The microINR Link Meter refers to the reader device.

The microINR Chips refer to the test strips, intended to be used exclusively with the microINR or microINR Link Meters.

The microINR System refers to reader (microINR or microINR Link Meters) and the strips (microINR Chips).

The microINR Link Meter can communicate through wireless Bluetooth technology to a software or smartphone application (App). Please, inform yourself on the correct use of this software or app.

Keep these instructions for use near the microINR Link Meter and refer to them if you have any questions about proper operation of the system.

At the end of these instructions you will find a glossary.

Self-testing or self-management of vitamin K antagonists may only be prescribed by a doctor.

Users of the microINR System (patients and healthcare professionals) must receive proper training before starting to use the system.

Safety Information

General Safety Warnings

Throughout these instructions for use you will find safety warnings and information on the correct use of the microINR System.

This warning symbol indicates a possibility of danger which could result in death, injury or harm to the patient or user if the procedures and instructions for use are not strictly followed.

This precaution symbol indicates the possibility of deteriorating or damaging the equipment and losing data, if the procedures and instructions for use are not strictly followed.

Important information regarding the correct use of the system that does not affect the safety of the patient or the integrity of the device is displayed over a blue background.

Infection Risk Control on Multi-Patient Test System

- Healthcare professionals must wear gloves during the entire process of the test.
- A separate lancet or lancing device must be used for each individual.
- Used Chips, lancets and gloves might be source of infection. Dispose of them in accordance with local regulations to prevent infections.
- Also, comply with your centre's internal hygiene and safety regulations.

There is a potential risk of infection. Healthcare professionals using the microINR system on multiple patients must take into account that all objects that come into contact with human blood are a possible source of infection. (See: Clinical and Laboratory Standards Institute: Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline, Fourth Edition CLSI document M29-A4, 2014).

Wireless Technology

The microINR Link Meter can communicate and transfer results to a PC or other electronic device (e.g., computer, tablet or smartphone) through a wireless Bluetooth connection.

Bluetooth is a form of radio frequency (RF) technology that

operates in the unlicensed industrial, scientific, and medical band at 2.5 GHz. The microINR Link Meter can only communicate with one compatible paired electronic device at each time; therefore, other devices with Bluetooth cannot access the information on the Meter unless the paired Meter is first unpaired and paired to a new electronic device.

In case of electromagnetic interference with another device, it is recommended to increase the distance between the Meter and that device. You could also turn off the interfering device. In addition, the Meter's Bluetooth wireless connection could be turned off.

The device generates, uses, and can radiate radio frequency energy.

Electromagnetic Compatibility and Electrical Safety

The microINR Link Meter system complies with electromagnetic compatibility requirements according to IEC 61326-1 and 61326-2-6 standards.

Do not use the microINR Link Meter near sources of intense electromagnetic radiation, as this could interfere with its correct operation.

The microINR Link Meter complies with electrical safety requirements according to standards IEC 61010-1 and IEC 61010-2-101.

The Meter must be disposed of as indicated in the WEEE Directive (Waste Electrical and Electronic Equipment 2012/19/EU).

1.3 ORAL ANTICOAGULANT THERAPY

Oral anticoagulant therapy is given to patients to prevent thromboembolic events such as venous thrombosis and pulmonary embolism or those linked to atrial fibrillation or artificial heart valves.

The treatment entails the need to monitor and adjust the doses periodically for each patient based on a blood test.

Depending on the reason for oral anticoagulant therapy a therapeutic range is defined for each patient meaning the value of the test should lie within that range.

Prothrombin Time and INR

The activity of oral anticoagulants is monitored by measuring the prothrombin time (PT) in seconds, which is the time it takes for a fibrin clot to form. The thromboplastin is used as a reagent to calculate the prothrombin time (PT). Depending on the nature of this reagent and the equipment used variations of the PT results are to be expected.

Therefore, in 1977, the World Health Organization (WHO) recommended a system standardisation method. Thus, prothrombin time values are converted into INR values, International Normalized Ratio, using the following equation:

$$INR = \left(\frac{PT}{MNPT} \right)^{ISI}$$

Where PT is the prothrombin time obtained in the test, MNPT is the mean normal prothrombin time and ISI the international sensitivity index corresponding to the thromboplastin.

The values of the MNPT and ISI parameters come from clinical calibration studies.

The pharmacological activity of vitamin K-antagonist oral anticoagulants can be modified by other drugs; therefore, you should only take the medicines prescribed by your physician.


Certain liver diseases, thyroid disorders and other diseases or conditions, as well as nutritional supplements, medicinal herbs or changes in diets can affect the therapeutic activity of oral anticoagulants and INR results.

1.4 MEASURING PRINCIPLE

The technology used by the microINR system is based on the microfluidics of the microINR Chip, which allows storing, dosing, moving and/or mixing small volumes of liquids to perform a chemical reaction.

The microINR Chips contain two channels, one for measurement and the other for control. An image of the Chip is shown below:

Each channel consists of a micro-reactor that contains the reagent and a microcapillary where the INR is determined. The reagent used in the measuring channel contains human recombinant thromboplastin and the reagent in the control channel contains recombinant thromboplastin and human coagulation factors to normalize the patient's blood.

The blood is applied to the Chip through the entry channel, separated into two channels and mixed with the reagents contained in each micro-reactor. The coagulation cascade is triggered instantly. When blood coagulates, its viscosity increases, which results into a change in blood flow behaviour.

The Meter captures the position of the sample by means of a Machine vision system and the position is transformed mathematically into speed and acceleration curves, from which an INR value is obtained.

Calibration

Each batch of microINR Chips has been calibrated against a reference batch of human recombinant thromboplastin traced to International Reference Thromboplastin of the World Health Organization¹.

These calibration values (ISI and MNPT) are encoded in the printed Datamatrix of each microINR Chip. Therefore, every test is automatically and individually calibrated eliminating any risk of human error.

(1). Expert committee of the WHO on biological normalization. Report forty-eight. Geneva, World Health Organization, 1999 (WHO technical report series No. 889)

2. microINR® SYSTEM

2.1 DESCRIPTION OF THE microINR® Link KIT

The microINR Link kit includes the following materials.
Make sure these items are included into the kit:

- Case
- microINR Link Meter
- Charger
- Plug adapter
- Instructions for use of the microINR Link Meter
- Quick guide
- Error guide
- Lancing device with instructions
(exclusive for Reference Kit KTD0001XX)
- Lancets (exclusive for Reference Kit KTD0001XX)

microINR Chips are sold separately.

Always carry the Meter inside its case.

2.2 PARTS OF THE microINR® Link METER

2.3 CHARGING THE microINR® Link METER

The Meter uses a lithium battery that can be recharged through the mini USB connection on the top of the Meter.

The recommended charging time is approximately 3 hours.

Charge the battery completely before using the Meter for the first time.

Do not open or manipulate the Meter.

Do not pierce or burn the battery.

Do not change the battery. The manufacturer will not warrant Meters that have been opened.

For battery replacement or meter repairs, the equipment must be sent to the manufacturer.

Use only the supplied charger and wall power adapter provided by the manufacturer or you may damage the Meter.

2.4 SETTING THE TIME & DATE

Before using the microINR System, check that the Meter's time and date are correct, since those are necessary to determine the Chips' expiration date.

Those parameters are set during the manufacturing process, so it probably will not be necessary to change them.

Time format: 24 hours. Date format: DD-MM-YYYY.

Follow the steps described below to change the Meter's time and date. Change only the digits needed to be set and confirm the rest digits:

- Press and hold the left and right buttons (E and M) at the same time for 10 seconds until the time field flashes.
- Press the left button (E) to set the hour.
- After selecting the correct hour, press the right button (M) and set the minutes.

- After selecting the minutes, press the right button (M) again and the date fields will start flashing.
- Use the left button (E) to select the correct day.
- When you reach the correct day, press the right button (M) to set the month. Use the left button (E) to select the correct month.
- When you reach the correct month, press the right button (M) to set the year. Use the left button (E) to select the correct year.
- After setting the month, time and date, press the right button (M) again to save your settings.

If the time and date setting is mistakenly initiated or a digit is unintentionally changed, the Meter will automatically exit that setting without saving any change after 10 seconds of inactivity.

2.5 PATIENT IDENTIFICATION (OPTIONAL)

microINR Link Meter can optionally record a result together with its Patient Identification (PID). To identify the patient, take the following steps before conducting a test:

- Hold the right button while inserting a Chip and release it afterwards. Upon releasing the button, "PID" and a 16-character alphanumeric field will be displayed.
- To enter the characters corresponding to the PID:
 - The first field will flash when selected.
 - To enter the first character, press the left button (E) until you reach the desired character.
 - Validate the character with the right button (M) and go to the next field.

- Repeat the previous steps until you reach the last field.
- After accepting the last field, the first field is automatically selected in case you want to correct any characters.
- To confirm the ID entered and go back, hold the right button (M) for 3 seconds.

While you are entering the PID, if one-minute elapses without pressing any button, an E01 message will be displayed.

When consulting the results on the Meter, the result of a test (INR or error message) is displayed along with its corresponding PID (if entered), date and time.

When entering the PID, do not include patient sensitive information (name, age, etc.) since it can be viewed in the Meter.

2.6 QUALITY CONTROL

Internal quality control

Meter performance is automatically checked when the system is turned on.

Integrated and independent on-board quality controls

Level 1 – Pre-test

- Chip integrity check
- Correct insertion check
- Automatic system calibration and rejection of expired Chips

Level 2 – Measuring channel

- Analytic verification performed on the measuring channel that identifies failures when processing the sample during the test, as well as proper pre-analytical treatment of the sample.

Level 3 – Control channel

- Control channel provides highly controlled clotting times. System reliability is assured when control clotting time lies within a pre-defined range.

PROFESSIONAL USE:

Liquid control:

The microINR system has a number of on-board quality control functions integrated into the Meter and the Chip and therefore there is no need to run quality control tests with liquid quality controls. However, iLine Microsystems has available an optional liquid control (plasma) for the microINR system. This control is provided to help meet the regulatory requirements applicable to your facility. To purchase, contact your local distributor.

3. CONDUCTING THE TEST

3.1 PREPARING THE NECESSARY MATERIALS

- microINR Link Meter
- microINR Chips CHB0001UU (not supplied)

- Fingersticking material
 - Professional use: disposable lancets (not supplied)
 - Self-testing: lancing device and lancets (exclusive Reference Kit KTD0001XX)
- Skin cleaning material (not supplied)


 Always use CE marked lancing devices and lancets.

3.2 MEASUREMENT PROCEDURE

Turning the Meter On

The Meter can be turned on:

- By inserting the Chip:
- By pressing any button:

Inserting the Chip

 Verify the expiry date and the storage conditions of the Chip before conducting the test.

- Open the pouch.

Hold the Chip by the yellow part so that the “microINR” logo can be read correctly. Insert the Chip into the slot and push it until it stops.


- Make sure the Chip has reached the end.

 If the pouch of the Chip is open, damaged or the film of the Chip is removed, please, dispose of the Chip and use a new one.

Modes of Use

The microINR system provides two application modes of the sample:


- Approaching the sample to the Meter.
- Approaching the Meter to the sample.

iLine Microsystems only recommends the option of approaching the Meter to the sample for users who are experienced in conducting INR assays with the microINR system.

Conducting the Test

- Once the Chip is inserted, the Meter will perform the quality controls mentioned in section 2.6 of these instructions for use. After conducting the aforementioned quality controls prior to the test, the “control” symbol will light up. If the first level of quality control is not passed, an error message will be displayed on the screen of the Meter.
- The Chip begins to flash and heats up until it reaches the proper temperature. Once this temperature is reached:
 - The device emits an audible signal (beep).
 - The drop symbol begins to flash on the display.
 - A count down appears (80s).
 - The Chip emits a steady light.
- Perform the fingerstick (see section 3.3 of these instructions for use).
- Make sure to obtain a spherical and properly sized drop, equivalent to a teardrop.
- Apply the drop of blood on the Chip immediately, by contact with the entry channel, without resting the finger on the Chip.
- The Meter will emit a beep when it detects the sample for the test and the drop symbol will stop flashing.

COUNTDOWN

- After the beep, gently move the finger away and wait until the INR results are displayed on the screen.

If you apply the sample and no sound is emitted, this means that there is not enough sample volume. Remove the Chip and repeat the test with a new Chip. Ensure that the size of the drop is enough and the entry channel is not blocked when placing it.

Do not reapply sample or try to add more blood to the channel.

Never perform the fingerstick before the start of the countdown.

Do not touch the Chip's entry channel with the finger while applying the sample.

Keep the Meter away from direct sunlight during the test.


Do not shake or drop the Meter. If the Meter is dropped or gets wet and the frequency of errors messages increases, contact your distributor.

Test Result and Assay End

- The measurement is performed and the Meter displays the result in INR units or an error message.


Error messages are displayed as a letter “E” followed by a number. If an error message is shown, follow the steps in the “Error Guide” section.

- Remove the Chip, holding it from both sides.

The used Chips, lancets and gloves might be source of infection.

For healthcare professionals: dispose the materials according to your institution's infection control policy and the appropriate local regulations. For self-testing patients: you can dispose all the materials in your regular waste. Discard the used lancets with care to prevent injuries.

Turning the Meter Off

There are two ways to turn the Meter off:

- The Meter turns off automatically after 3 minutes of inactivity.
- Press the left button (EXIT) to turn the Meter off, holding it for 3 or 4 seconds.

The Meter cannot be turned off while connected to the power supply.

3.3 OBTAINING AND APPLYING THE CAPILLARY BLOOD SAMPLE

The steps to **obtain and apply a capillary blood sample correctly** are detailed below:

- Read the instructions for use of the fingerstick device or lancet.

Before lancing the finger, it is convenient to warm hands. There are several techniques that can be used for that purpose such as washing hands with warm water, keeping hands below the waist, massaging the fingertip softly...

The fingerstick site must be clean, completely dry and free of contaminants. Washing the hands with warm soapy water is recommended. You may also use alcohol to clean the fingerstick area. Always dry the area thoroughly to remove any traces of substances that might interfere with the result. Always use a new, clean and dry gauze.


Any alcohol (disinfectants, shaving creams, etc.), lotions or sweat on the fingerstick area or the blood sample may cause incorrect results.

- You can use any finger for the fingerstick. The recommended site is the one shown on the following image.
- Once the Meter is ready for testing and the 80-second countdown appears on the screen, place the lancet firmly against the finger and press the button. Press the finger

gently from the base to help form a drop of blood. Do not press the fingerstick site or let the drop of blood spread on the finger.

- Before placing the sample on the Chip, make sure to obtain a spherical and properly sized drop (equivalent to a teardrop), large enough to leave a small amount of blood (remnant) at the entry channel.


Sampling technique can affect the result of the test. Do not squeeze or "milk" the fingerstick area as this can alter the coagulation process.

When using the Meter in a multi-patient site uses only auto-disabling single use lancets.

- Apply the drop to the Chip immediately, in contact with the entry channel.


Samples must be applied immediately after collection, since blood clotting does naturally occur upon fingerstick. Avoid contact between the Chip and the finger in order not to obstruct the entry channel and to allow for uninterrupted blood absorption. Only the drop of blood must make contact with the Chip.


Apply the sample on a single attempt. Never add more blood to the Chip.

- After the beep tone, gently remove the finger, leaving a small amount of blood (remnant) at the entry channel as shown in the picture.
- If the test must be repeated, perform the fingerstick on a different finger, using a new Chip and lancet.

3.4 INTERPRETING THE RESULTS

The results are shown as International Normalized Ratio (INR) units. The microINR System's results range between 0.8 and 8.0.

If you obtain a result out of the measuring range, the display will show ↓ 0.8 (below 0.8) or ↑ 8.0 (above 8.0).

If an error message is displayed, see the “Error guide” section and follow the instructions.

If an unexpected result is obtained, repeat the test making sure that the indications described in these instructions for use are strictly followed. If an unexpected result is obtained again, contact your physician and/or distributor.

Results are unexpected when they lie outside the therapeutic range or do not match the patient’s symptoms: haemorrhages, bruises, etc.

3.5 LIMITATIONS OF USE

- Once the Chip is taken out of its original pouch, it should be used within the following 6 hours.
- Do not use to measure or monitor the anticoagulation status of patients under treatment with direct oral anticoagulation treatments (non-vitamin-k antagonist drugs).
- The performance of the microINR System has not been demonstrated on blood samples with hematocrit values outside the range of 25% to 55%. Hematocrit out of this range may affect test results.
- The device is highly sensitive to vitamin K dependent coagulation factor deficiencies.

For more information on interferences of the microINR System with other drugs and diseases and the microINR System performance, refer to the instructions for use of the microINR Chip.

4. MEMORY

The microINR Link Meter can store up to 199 results of patients, liquid quality controls and error messages. Each result is stored with the date and time of the test. When conducting a test, if there is no free storage space, the oldest result will be automatically deleted to store the new one.

To check the results:

- Press the right button (M). The result of the last test conducted will be displayed with its date and time and the ID of the corresponding patient (if it was entered).
- Press again to display the next result, corresponding to the second-to-last test and so on.
- Press the left button (E) to return to the initial screen. If you enter a Chip while you are checking the memory, a new test will begin normally.

5. BLUETOOTH® FUNCTIONALITY

The microINR Link Meter can communicate and transfer results to a PC or other electronic devices through wireless Bluetooth connection.

Transferring results from microINR Link Meter is useful for data storing in DMS (Data Management Software) or in suitable smartphone/tablet Apps.

Please, contact your local distributor for further information about DMS or Apps compatible with the microINR Link Meter. microINR Link Meter can just be paired with one device at a time.

A description of all the symbols related to the microINR Link Meter Bluetooth performance can be found at the end of this instructions of use or at the Quick Guide included in the microINR Link Kit.

In case the “Bluetooth failure” symbol displays on the Meter’s screen, contact with your local distributor.

5.1 BLUETOOTH® CONFIGURABLE FEATURES

The microINR Link Meter wireless communication via Bluetooth can be configurable to suit the user’s needs. Contact your local distributor if you want to change the default microINR Link Meter configuration. The configurable features are:

- The microINR Link Meter Bluetooth functionality can be turn on or off.
- Pairing modes:
 - Free mode: no password required.
 - Password mode: a pairing code will be required to connect both devices.
 - NFC Mode: allows pairing devices by approximation.
- If, after having been paired with a device, the microINR Link Meter is paired with another one, it is possible to configure the management of the results performed on the microINR Link Meter while paired with the previous device.


5.2 microINR® LINK METER BLUETOOTH® MENU

The microINR Link Meter Bluetooth functionality has been designed to be as simple and automatic as possible. For certain processes, it will be necessary to enter the microINR Link Meter Bluetooth menu.

To enter the microINR Link Meter Bluetooth menu, ensure that the Meter is turned off and no cable is connected to the USB port. Turn the Meter on by pressing the left and right buttons at the same time. Hold both buttons (approx. 4 seconds) until the Meter shows the menu’s first option.

The menu allows to perform the following four activities:

- 1. Pairing:** to pair the microINR Link Meter, see the section 5.3.

2. Unpairing: to unpair the microINR Link Meter, see the section 5.4.

3. Send all results: to send all the result stored at the microINR Link Meter, see the section 5.6.

4. Exit.

To move through the different menu options, press the right button. To execute the desired option, press the left button. You can exit the Bluetooth menu at anytime by inserting a microINR Chip or after 20 seconds of inactivity the Meter will skip the menu automatically.

5.3 DEVICE PAIRING

• The microINR Link Meter can be paired with a compatible device automatically or through the Bluetooth Menu:

- When the microINR Link Meter is turned on for the first time by pressing any button, it, automatically and for a minute, will try to pair with a compatible device.


- The pairing process can also be initiated through the Bluetooth menu (see section 5.2).

• Different steps must be followed depending on the pairing mode set at the microINR Link Meter:

- Free mode: just follow the instructions on the compatible software or smartphone/tablet App.

- Password mode: the microINR Link Meter shows a 6-digit password. Type the password on the compatible device and follow the software or App instructions.

- NFC mode: approximate the microINR Link Meter top back side with the compatible device so that they can get in contact. Follow the software or App instructions.

- When the microINR Link Meter is paired the “*Process Correct*” symbol appears on the Meter’s screen. The microINR Link Meter can only be paired with one compatible device at a time. Unpair the microINR Link Meter before pairing with another device (see section 5.4). If a new pairing process is attempted while the microINR Link Meter is already paired, the “*Process Incorrect*” symbol will be displayed on the Meter’s screen for 4 seconds.


Once the pairing process starts, the microINR Link Meter is accessible for 60 seconds to pair with a compatible device. If the microINR Link Meter does not connect with another device, the “*Process Incorrect*” symbol will appear. The “*Unpaired Device*” symbol will appear on the lower right side of the screen.

The Pairing process can be cancelled at any time by pressing the left button or inserting a Chip.

5.4 DEVICE UNPAIRING

- Enter the Bluetooth Menu to unpair the microINR Link Meter. Ensure that the Meter is turned off and no cable is connected to the USB port. Turn the Meter on by pressing the left and right button at the same time for 4 seconds.
- Press the right button to reach the unpairing option and then press the left button.
- The unpairing process starts automatically. After a few seconds the “*Process Correct*” symbol is displayed on the Meter’s screen for 4 seconds. The “*Unpaired Device*” symbol appears on the lower right side of the screen.

To skip the Bluetooth Menu, press the right button until you get to the “*Exit*” option (see section 5.2). You can also insert a microINR Chip or after 20 seconds of inactivity the Meter will skip the Menu automatically.

Once the microINR Link Meter is unpaired, it is advisable unpairing also the electronic device that was previously connected to the Meter.

5.5 TRANSFERRING DATA

Once the microINR Link Meter is linked with a compatible device, it automatically sends the result or the error message obtained after each test. During the data transferring process, the “Bluetooth” symbol appears on the lower right side of the screen and it blinks while sending the result.

If successfully sent, the “Result sent” symbol will be displayed on the Meter’s screen for 4 seconds.

If the linked device is not available, the microINR Link Meter will store the result and this will remain pending to be sent. The “Not Sent result” symbol will be display.

In case of having results pending to be sent, the “Stored-pending result” symbol appears on the lower right side of the display. The microINR Link Meter will try to send those pending results automatically every time the Meter is turned on and/or after performing a new test.

5.6 SENDING ALL RESULTS

It is also possible to send all results (up to 199) stored in the microINR Link Meter memory.

- To send all stored results, enter the Bluetooth Menu. Ensure the Meter is turned off and no cable is connected to the USB port. Turn the Meter on by pressing the left and right button at the same time for 4 seconds. Press the right button until you reach the “Send all results” option and then press the left button.
- The Sending process starts automatically. During the sending process the following symbol appears. When the result’s transfer is finished the “Results sent” symbol appears on the Meter’s screen.
- To skip the Bluetooth Menu, press the right button until you get to the “Exit” option (see section 5.2). You can also insert a microINR Chip or after 20 seconds of inactivity the Meter will skip the Menu automatically.

If the microINR Link Meter is paired with a new device, the previously stored results can be protected, and their transfer avoided. Please contact your local distributor to help you with the microINR Link Meter configuration.

6. CLEANING AND DISINFECTING THE microINR® Link METER

Cleaning and disinfection of the microINR Link Meter is essential to ensure proper microINR System operation and to prevent blood-borne transmission of pathogens in multi-patient tests.

Clean the Meter when it is visibly dirty and before disinfecting. Use a new pair of gloves every time you clean and disinfect the Meter.

Before cleaning or disinfecting the Meter, turn it off and make sure the cables are unplugged.

To clean the Meter:

- Clean the Meter with a clean gauze or wipe moistened with isopropyl alcohol 70% until there is no visible dirt.
- Ensure there are no remaining fibers or lint on any part of the Meter, especially on the Chip insertion area and the USB port.

To disinfect the Meter:

- Disinfect all the parts of the Meter with a clean gauze or wipe moistened with isopropyl alcohol 70%.
- Wait one minute for the alcohol to act.
- Thoroughly dry the Meter with a dry and clean lint-free cloth or gauze.
- Afterwards, wait 15 minutes for the Meter to dry and ensure it is completely dry before conducting a new test.

- Ensure there are no remaining fibers or lint on any part of the Meter, especially on the Chip insertion area and the USB port.
- Discard the used gauze, wipes and gloves.

Do not clean or disinfect the Meter while conducting a test.

Do not use aerosols or any cleaning or disinfecting agents other than a clean gauze or wipe moistened with isopropyl alcohol 70%.

Verify that the gauze or wipe is just moist, not soaked.

Do not spray fluids on the Meter or submerge the Meter.

Make sure that no fluids enter the Meter or the Chip insertion area.

The Chip insertion area must always be clean and dry before conducting a test. When inserting a Chip, remains of blood or alcohol can contaminate the sample.

Do not handle the Chips with alcohol-contaminated hands/gloves.

Comply with all recommendations regarding cleaning and disinfection of the Meter.

Not doing so could cause incorrect results.

7. ERROR GUIDE

Error	Probable Cause	Possible Solution
Messages when preparing to test		
E01	The Datamatrix could not be read.	Insert the same Chip again, ensuring correct insertion. If the problem persists, repeat the test with a new Chip.
E02	Expired Chip.	Verify the date of the Meter. If the date is not correct, enter the current date and insert the same Chip again. If the date is correct, repeat the test with a new batch of Chips after verifying their expiry date.
E03	The 80 second interval for applying the sample has been exceeded.	If the sample has not yet been applied, repeat the test with the same Chip.
E04	Chip inserted backwards.	Rotate the Chip and repeat the test.
E11	Faulty or incorrectly inserted Chip.	Insert the same Chip again, ensuring complete insertion. If the problem persists, repeat the test with a new Chip.
Messages related to the test		
E05/ E09	Inadequate coagulation of the sample during the test.	Repeat the test with a new Chip. Strictly follow instructions on obtaining and applying the sample. If the problem persists, repeat the test with a new box or batch of Chips.
E10	The INR value of the control channel is outside the defined range.	Repeat the test with a new Chip. If the problem persists, repeat the test with a new box or batch of Chips.
E14/ E15/ E17	Error while processing the sample during the test.	Repeat the test with a new Chip. Strictly follow instructions for conducting the test. If the problem persists, repeat the test with a new box or batch of Chips.
E16	Inadequate coagulation of the sample during the test. ATTENTION: Possible sample with abnormally high clotting times.	Repeat the test with a new Chip. If error E16 is displayed again, use a different measurement method. This error occurs mainly in patients with high clotting times.
E18	Inadequate sample handling or hematocrit outside the defined range.	Repeat the test with a new Chip. Strictly follow instructions on obtaining and applying the sample. If error E18 is displayed again, use a different measurement method. This error occurs mainly in patients with a hematocrit outside the defined range for the microINR system (25%-55%).
Other messages		
E06	Failure while checking the electronic components of the Meter.	If the problem persists, contact your local distributor.
E07	Temperature below the defined range.	Repeat the test in a warmer location.
E08	Low battery.	Charge the device with the charger supplied by the manufacturer.
E12	Temperature above the defined range.	Repeat the test in a cooler location.

8. ADDITIONAL INFORMATION

8.1 SPECIFICATIONS

- Dimensions of the Meter: 119 x 65 x 35 mm.
- Weight: 225±3 gr. (Battery included).
- Screen: LCD 45 x 45 mm.
- Memory: 199 results / error messages with their date and time.
- *Bluetooth*[®] low energy 5.0 backwards compatible with 4.2 version*.
- Power supply:
 - Battery: Lithium 2400 mAh/2800 mAh; 3.7 V.
Consumption: 1 A.
 - Power supply: Only plug the charger to power supplies with the following characteristics (Input): 100-240 V, 50-60 Hz,
Consumption: 0.6 A.
 - Power supply (Output): 5 V dc, 1200 mA via a mini USB connection.
- Battery life: **approximately 70 tests.
- Operation conditions:
 - Temperature: 15 °C - 35 °C.
 - Maximum relative humidity: 80%.
- Meter storage temperature: -20 °C to 50 °C.

- Measurement range: 0.8 – 8.0 INR.
- Sample volume: at least 3 µL.
- Data transfer via a mini USB connection*.
- * For data transfer via USB or Bluetooth connection, please, contact your local distributor.
- ** Test conducted at 22 °C with a 10-minute period between tests.

8.2 WARRANTY

iLine Microsystems warranties to the original buyer that the microINR system is free of material and manufacture defects for two years after the purchase date.

This warranty does not cover any component damaged due to inadequate storage in environmental conditions outside the defined range, accidents or modifications, incorrect use or handling and misuse. The buyer must deliver a written warranty complaint to the manufacturer within the corresponding warranty period.


8.3 TECHNICAL SERVICE


If a problem persists after performing the actions stated in the table of errors or if you want additional information, you can contact your local distributor.

Note that any repair or modification to the microINR System can only be performed by staff authorized by iLine Microsystems.


8.4 SYMBOLS

	"Manufacturer"
	"Serial number"
	"Catalogue number"
	CE "CE Marking"
0123	0123 Notified Body Identification Number <i>(only applies for the Selftesting Use certification)</i>
	"Electronic waste selective collection"
	"Consult the instructions for use"
	"in vitro diagnostic medical device"
	"Direct current"
	"Temperature limit"
	"Batch code / Lot Number"
	"Biological hazard"
	"Warning"
	"Precaution"

Menu:
Pairing


Menu:
Unpairing

Menu:
Send all
results

Menu:
Exit

Pairing:
Free Mode

Pairing:
Password

NFC

Pairing:
NFC

Unpaired
Device

Sending

Result/s
sent

Result/s
not sent

Pending
result/s

Bluetooth
Chip on

Bluetooth
failure

Process
Correct
Incorrect

8.5 GLOSSARY

Bluetooth: is a wireless technology standard for exchanging data between fixed and mobile devices over short distances using short-wavelength UHF radio waves in the industrial, scientific and medical radio bands, from 2.400 to 2.485 GHz, and building personal area networks (PANs).

Capillary blood: blood from the smallest blood vessels in the body, usually obtained by puncturing a fingertip.

Capillary fingerstick: small puncture on a finger to obtain capillary blood.

Chip: disposable element inserted into the Meter. Serves to enter the sample for the INR test.

Control channel: channel used to measure the normalised coagulation time in order to detect degradation of the reagents.

Entry channel: slot on the bottom of the Chip that receives the blood.

International Normalized Ratio [INR]: standardised prothrombin time measurement system that accounts for the different sensitivity of the thromboplastins used in different methods. The INR results from different prothrombin time measurement systems can be compared to each other.

Lancet: piercing tool used to make a small cut or puncture to collect a small drop of blood, which will be used for the INR test.

Microcapillary: site on which the INR test takes place.

Microfluidics: technology for storing, dosing, transferring and/

or mixing small volumes of fluid to cause a chemical reaction.

MicroINR Link Meter: electronic device that serves to conduct INR tests.

Microreactor: area of the Chip meant to store the reagents.

Mini-USB connector: connector on the top front of the microINR Link Meter.

Oral anticoagulant therapy: orally administered treatment that inhibits or interferes with the coagulation of the blood.

Plasma: liquid part of the blood.

Prothrombin time (PT): coagulation analysis performed to analyse the extrinsic coagulation pathway.

Quality control: tests to prove that a system performs correctly providing reliable results.

Reagent: substance used to cause a chemical reaction in order to measure a substance or process (such as the INR test).

Remnant: small amount of excess blood that remains on the entry channel of the Chip.

Therapeutic range: safe range of INR values. The physician determines a specific therapeutic range for each patient.

Thromboplastin: A substance used by blood platelets and combined with calcium that converts prothrombin (protein) into thrombin (enzyme) as part of the clotting cascade.

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